

## WHAT KITS ARE RIGHT FOR YOU?



### ALVEOSAMPLER™

The AlveoSampler system is a disposable collection system used to collect alveolar air in a standard syringe for immediate analysis, and is the most commonly used.

Single-patient use of this device eliminates the danger of inter-patient cross-infection, and will save time and money related to the cost of cleaning and sterilizing reusable collection components.

This system can only be used in-house on any BreathTracker™ or MicroLyzer.

Each kit will contain all the necessary supplies and substrate desired for collection.



### GASAMPLER™

The GaSampler system is used for the sampling and storing of alveolar air for immediate or later analysis. This system is designed for only in-house testing on any BreathTracker or MicroLyzer.

This system can be used by non-technical personnel, or even by the patient without supervision after an explanation of the procedure.

This system uses multiple parts that do not come together in a kit; all parts are sold as disposable except the multi-patient collection bag which can be used multiple times over multiple patient collections.

These collection bags can be custom made to any size breath sample desired\*.



### KIDSAMPLER™

The KidSampler system is similar to the GaSampler, except it can utilize a mouthpiece or a face mask for pediatric collection. This system is designed only for in-house testing on any BreathTracker or MicroLyzer.

This system uses multiple parts that do not come together in a kit; all parts are sold as disposable except the multi-patient collection bag which can be used multiple times over multiple patient collections.

These collection bags can be custom made to any size breath sample desired.

Other pediatric collection systems are available for babies and neonatal patients.



### EASYSAMPLER™

The patented EasySampler allows patients to collect samples off-site unsupervised for later analysis or can be used in-house as well, but not recommended for patients who cannot follow verbal directions.

Kits are provided with easy instructions and all necessary supplies for the patient to perform the test. Each kit can be customized with your own literature and mailing labels\*.

This system can only be used with the purchase of a BreathTracker SC or H2 Plus analyzer and SamplXtractor™ system.

Samples can be stored for up to 14 days with no loss in concentration, and are designed for easy and safe transport.