



FREE UP CONGESTION IN YOUR OFFICE, SEND PATIENTS HOME WITH TEST KITS!



The SamplXtractor™ Model SX-2 Plus is a stand-alone accessory designed for use with the BreathTracker™ and MicroLyzer™ breath-testing instrumentation.

Simplify the extraction of samples from vacuum tubes that are used with the patented EasySampler™ system for reliable injection into your breath-testing instrumentation.

By utilizing the SamplXtractor system, you now have the ability to use EasySampler™ kits to send patients off-site with test kits or transport samples between sites. This is especially useful if you cannot analyze samples immediately or cannot afford to congest your office with multiple patients testing over the length of time needed for sample collection.

Use of this system is extremely simple and ensures the best extraction of patient samples from QuinTron evacuated tubes. Minimal maintenance requirements and no necessity for calibration make this accessory extremely beneficial and easy to use and maintain for nurses/staff or technicians of any type.

Breath Collection Products that are for use with the SamplXtractor Model SX-2 Plus are:

EasySampler System

The patented* EasySampler allows patients to collect samples off-site or in-house unsupervised for later analysis and is only recommended for patients that can follow simple verbal and written instructions.

Kits are provided with easy instructions and all necessary supplies for the patient to perform the test and are designed for easy and safe transport.

Each kit can be customized with your own literature, protocol, substrate, and mailing labels.†

Samples can be stored for up to 14 days with minimal to no loss in sample integrity.

Catalog Number: QT01200-EX

Specifications:

Width: 6.5"
Height: 5.5"
Depth: 8"
Weight: 3.75 lbs (1.70 kg)
Power: 24 Volts D/C

TAKE HOME KITS AVAILABLE:

- Lactose
- Bacterial Overgrowth with Lactulose
- Bacterial Overgrowth with Glucose
- Fructose
- Sucrose
- d-Xylose
- Custom Kit development available



* US Patent # 5,467,776
EasySampler Test (Substrate Not Shown)



3712 West Pierce Street, Milwaukee, WI, USA
www.QuinTron-USA.com
Phone: (800)542-4448 / (414)645-4222

† Customized test kits are available upon request and have development fees not included in the cost of standard kits.